



Milkin' Mamas Breast Milk Donations

Breast Milk Donation Check List:

- Make a personal commitment to completing the donor application process and making a breast milk donation.*
- Complete the Preliminary Questions for potential donors, either via email or over the phone.*
- Make sure you understand the donor application process and have had all of your questions answered.*
- Complete and submit the Prolacta Bioscience on-line Donor Application.*
- Respond to any emails asking for additional information in order to approve your application (if applicable).*
- Print, sign and return the Medical Confirmation Forms emailed to you.*
- Upon receipt of the Donor Test Kit (in this order):*
 - *Email your freezer temperature reading, using thermometer provided.*
 - *Email the last 5 digits printed underneath the bar codes on the DNA swabs and blood vials.*
 - *Perform the DNA collection using the swabs provided and return in white prepaid FedEx envelope.*
 - *After an order has been placed, schedule your blood draw. Mon-Thurs in the morning or early afternoon.*
 - *Place the small gel packs provided with the test kit in your freezer.*
 - *Examiner will perform your blood draw with the supplies provided and ship your blood sample to the lab.*
- Once Qualified as a donor, reply to our email and provide the following information:*
 - *How many ounces and bags of milk you have ready to donate.*
 - *Are you making a one-time donation or do you plan on being an on-going donor.*
- Upon receipt of the shipping cooler:*
 - *Place the gel bricks provided in your freezer for at least 48 hours.*
 - *Ship your milk donation preferably Monday through Wednesday, NEVER on a Friday.*
 - *Pack your milk donation according to the instructions provided.*
 - *Call FedEx (1-800-GO-FEDEX) to pick up your overnight delivery prepaid package.*
 - *Email us (donate@milkinmamas.com) with the total number of ounces and bags you were able to share.*
- Congratulate yourself on making a difference in the lives of premature babies and giving hope to their families.*
- If you are an on-going donor, request additional shipping coolers and milk storage bags as needed.*

As always, never hesitate to contact us if you have any questions or concerns. We are here to help!

www.MilkinMamas.com - donate@milkinmamas.com - 562-421-6969